

ri-thermo® fastPRObe

Riester Fast Reading Thermometer

- › Fast and reliable spot-check measurements
- › Riester, colour coded probe to distinguish between oral/axillar (blue) and rectal (red), used with probe covers to reduce risk of cross infection
- › 1000 + measurements with one set of batteries
- › Switchable from Celsius to Fahrenheit
- › Automatic On/Off



Measurement progress bar



Convenient storage for 20 probe covers

Backlight illuminated and easy to read LCD screen

Location button for optimized measurements

Pulse timer

Easy-to-use recall button to review the last measurement

Professional electronic thermometer for quick and accurate body temperature measurements. Designed for professional use in the daily routine. Support for oral, rectal and axillary measurements.

Medical equipment:	Electronic thermometer
Electrical protection:	Device with internal voltage supply
Electrical degree of protection:	Type B
Operating mode:	Continuous manual operation
Dimensions / Weight:	D: 205 x 70 x 70 mm / W: 400 g with batteries and probe cover box
Battery	4 x AA, LR6 Alkaline 1.5 V or 4 x AA Mignon NiMH 1.2 V
Temperature measurement range:	34°C to 43°C (93.2°F to 109°F)
Ambient operating temperature range:	10°C to 40°C (50°F to 104°F), at 10 % - 85 % non-condensing
Storage and transport conditions:	-20°C to 55°C (-4°F to 109°F), at 10 % - 95 % non-condensing
Air pressure	700 to 1050 hPa
Laboratory accuracy water bath accuracy (36° to 42.0°C):	Direct measuring mode ± 0.1°C (± 0.2°F), fast measuring mode ± 0.3°C (± 0.5°F)
Direct / Fast measuring mode accuracy (all measuring points)	D: ± 0.1°C (± 0.2°F) / F: ± 0.3°C (± 0.5°F)
Typical measurement times (after position of the probe)	Direct measuring mode oral mode: 60 sec., axillary mode: 60 sec., rectal mode: 60 sec. Fast measuring mode oral mode: 12 sec., axillary mode: 12 sec., rectal mode: 12 sec.

Art.No. 1830	ri-thermo® fastPRObe Thermometer blue (oral, axillary)
Art.No. 1831	ri-thermo® fastPRObe Thermometer red (rectal)
Art.No. 13335	Thermometer probe - probe insert blue (oral, axillary)
Art.No. 13336	Thermometer probe - probe insert red (rectal)
Art.No. 13334	Probe cover package containing 20 pieces
Art.No. 13338	Probe cover package containing 25 x 1 package of 20 pieces
Art.No. 12955	Calibration kit for ri-thermo® fastPRObe